DK50 DE

- EN User manual
- DE Benutzerhandbuch
- RU Руководство пользователя
- SK Návod na použitie

MEDICAL COMPRESSOR
MEDIZINALKOMPRESSOR
МЕДИЦИНСКИЙ КОМПРЕССОР
MEDICINÁLNY KOMPRESOR

DK50 DE

EKOM spol. s r. o.
Priemyselná 5031/18
SK-921 01 Piešťany
Slovak Republic
tel.: +421 33 7967255
fax: +421 33 7967223

www.ekom.sk
email: ekom@ekom.sk

DATE OF LAST REVISION
DATUM DER LETZTEN ÜBERARBEITUNG
ДАТА ПОСЛЕДНЕГО ПЕРЕСМОТРА
DÁTUM POSLEDNEJ REVÍZIE

04/2026

NP-DK50 DE-LF-17_04-
2026
112000082-000

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CONTENTS

1.	GENERAL INFORMATION	6
	INTENDED USE	6
	CONTRAINDICATIONS AND SIDE-EFFECTS	6
	OPERATOR'S RESPONSIBILITY FOR PATIENT SAFETY	6
	MARKINGS.....	6
	WARNINGS	6
	General safety warnings	7
	Electrical system safety warnings.....	7
	WARNING NOTICES AND SYMBOLS	8
	USE.....	9
	STORAGE AND TRANSPORT	9
2.	EQUIPMENT DESCRIPTION	10
3.	TECHNICAL DATA	12
	FAD efficiency correction for differences in elevation	13
4.	OPERATION.....	14
	INSTALLATION AND FIRST OPERATION.....	14
	Removal of transport stabilizers	15
	Compressed air connection.....	15
	Electrical connection.....	15
	First operation.....	16
	Adjusting the output air pressure	16
	OPERATION.....	16
	Switching the compressor on.....	16
	Running the compressor	16
	Cleaning the compressor	17
5.	MAINTENANCE	18
	REPAIRS AND SERVICE.....	18
	Cover removal	18
	MAINTENANCE SCHEDULE	19
	Checking tightness of joints and inspecting the equipment	19
	Replacing filter elements	20
	Changing suction filter elements	21
	Securing the equipment before transport	21
	SHUTDOWN	21
	EQUIPMENT DISPOSAL.....	21
6.	TROUBLESHOOTING.....	22
7.	SPARE PARTS	23
8.	ELECTRIC AND PNEUMATIC DIAGRAMS.....	24
	WIRING DIAGRAM.....	24
	PNEUMATIC DIAGRAM.....	25
9.	ELECTROMAGNETIC COMPATIBILITY DECLARATION	26

GENERAL INFORMATION

1. GENERAL INFORMATION

INTENDED USE

The EKOM DK50 DE is a medical air compressor that supplies clean, oil-free compressed air for use with medical ventilators.

CONTRAINDICATIONS AND SIDE-EFFECTS

There are no contraindications or side-effects known.

OPERATOR'S RESPONSIBILITY FOR PATIENT SAFETY

Instructions for use is an integral part of the equipment and must be kept with the compressor. Careful review of this manual will provide information necessary for correct operation of the equipment.

Rx only

US Federal law restricts the sale of this device by or on the order of a physician.

MARKINGS

Products marked with the CE mark of compliance meet the safety requirements of the European Union (93/42/EEC).

WARNINGS

- The safety of operating personnel and trouble-free operation of the equipment are ensured only if original parts are used. Only optional equipment and spare parts mentioned in the technical documentation or expressly approved by the manufacturer can be used.
- If any other optional equipment or consumable materials are used, the manufacturer cannot be held responsible for the safe operation and functionality of the equipment.
- The warranty does not cover damages resulting from the use of optional equipment or consumable materials other than those recommended by the manufacturer.
- The manufacturer assumes responsibility for the safety, reliability and function of the equipment only if:
 - Installation, calibration, amendments, extensions and repairs are performed by the manufacturer, one of its representatives or a service provider authorized by the manufacturer
 - The equipment is used in accordance with the Installation, Operation and Maintenance Manual
- The Installation, Operation and Maintenance Manual accurately describes the design of the compressor and its compliance with safety and technical standards. The manufacturer reserves all rights to its wiring diagrams, procedures and names.

- This user manual is the original instructions. Translation is performed in accordance with the best available knowledge.

General safety warnings

The equipment is designed to operate safely when used correctly. Please note the following safety measures to avoid injury or damage.

- Equipment operation must comply with all local codes and regulations.
- Original packaging should be kept for the possible return of the unit. Only original packaging ensures optimal protection of the equipment during transport. If it is necessary to return the equipment during the warranty period, the manufacturer is not liable for damages caused by incorrect packaging.
- The user must immediately notify the supplier if any problem occurs during the use of the equipment.
- This product is not intended for use in areas where there is a risk of an explosion. Do not operate the compressor in the presence of flammable anesthetics.
- Never feed oxygen or nitrous oxide into the compressor. Compressor components are not approved for oxygen or nitrous oxide use.
- The compressor must not be used for supplying air to the medical air central distribution pipeline system.

Electrical system safety warnings

- The equipment must be connected to ground. In order to assure proper grounding, connect the compressor to a receptacle marked "hospital grade."
- Before the compressor is plugged in, make sure that the voltage and frequency of the mains specified on the equipment are the same as the power mains.
- Before operating, check for possible damage to the equipment and any connections. Damaged pneumatic and electrical lines must be replaced immediately.
- If a technical failure occurs, immediately disconnect the equipment from the mains (pull out the main power plug).
- During repairs and maintenance, ensure that:
 - The main power plug is removed from the power socket
 - Compressed air lines are disconnected
 - All pressure has been released from the air tank
- Only a qualified technician can install this equipment.

GENERAL INFORMATION

WARNING NOTICES AND SYMBOLS

The following symbols are used for important information in the Installation, Operation and Maintenance Manual and on packaging and the product:

Attention, see instructions for use

Caution, risk of electric shock

Refer to instruction manual

CE mark of compliance

Medical device


Serial number

Article number

Unique Device Identifier

Swiss Authorised Representative

Swiss Importer

Compressor is controlled automatically and may start without warning

Caution, hot surface

Handling mark on package – FRAGILE

Handling mark on package – THIS SIDE UP

Handling mark on package – KEEP DRY

Handling mark on package – TEMPERATURE LIMITATIONS

Handling mark on package – LIMITED STACKING

Mark on package – RECYCLABLE MATERIAL

Ground connection

Equipotentiality

Fuse

Condensate drain

Never dispose of the battery with common household waste.

USE

- The equipment can be installed and operated only in a dry, ventilated and dust-free area. Climatic conditions for operation - see Technical data.
- The compressor must stand on a flat and stable base.
- The compressor must not be exposed to rain. The equipment must not be used in humid or wet environments. Never use the compressor in the presence of flammable liquids or gases.
- Before connecting the compressor to respiration equipment, make sure that it meets the requirements of the respiration equipment. Refer to the Technical data for this purpose.
- Any use other than the compressor's intended use is not considered to be safe. The manufacturer is not responsible for any damages that result if the compressor is used for any other purpose. Risk is exclusively assumed by the operator/user.

STORAGE AND TRANSPORT

The compressor is shipped from the factory in transport packaging with the pump stabilized, protecting it from damage during transport.

For transport, always use the original packaging and secure the compressor in the upright position.

Protect the compressor from humidity, contamination and extreme temperatures during transport and storage. A compressor in its original packaging should be stored in a warm, dry and dust-free area.

Keep the packaging material, if possible. If not, dispose of the packaging material in an environmentally-friendly way. Cardboard can be recycled.

Prior to transport it is necessary to secure the motor inside the compressor (Chapter 5).


EQUIPMENT DESCRIPTION

2. EQUIPMENT DESCRIPTION

1. OUT - compressed air outlet
2. Pressure regulator
3. Output air pressure gauge (with accuracy of $\pm 5\%$)
4. Hour meter
5. Switch, power connector and fuses
6. Power cord anchor
7. Equipotential connector
8. Compressor
9. Suction filter
10. Fan
11. Safety valve
12. Filter (5 μm) with automatic condensate drain
13. Supporting base
14. Cooler
15. Regulator safety cover

The equipment uses an oil-free piston compressor (8) driven by a low-maintenance electric motor. Compressed air is cooled in the cooler (14) where condensed water is separated into a separate tank (12). Compressor working pressure is set using the safety valve (11). Constant output pressure is maintained by the pressure regulator (2).

EQUIPMENT DESCRIPTION

TECHNICAL DATA

3. TECHNICAL DATA

TYPE		DK50 DE
VERSION		low flow
Output flow at pressure 3.5 bar (51 psig)	L/min	20
Voltage / Frequency / Rated current	V/Hz/A	230/50 / 2.3 230/60 / 2.3 110-120/60 / 4.4
Air filtration	µm	5
Pressure dew point at 3 bar, 20°C, 50% Relative air humidity		5°C (9°F) below the ambient temperature
Outlet connection		DISS 1160-A (3/4"-16 UNF)
Sound level ^{a)}	dB(A)	48 49 / 60 Hz
Mode of operation		Continuous - S1
Separation of condensed water		Automatic
Output pressure		3.0 bar (43 psig) Adjustable to max. 3.5 bar (51 psig)
Pressure range		5 bar (72.5 psig)
Adjustment of output air pressure		Pressure regulator
Dimensions of compressor	w x d x h	395 x 270 x 400 mm (10.5 x 15.5 x 16 in)
Dimensions of packaging	w x d x h	460 x 350 x 500 mm (14 x 18 x 19.5 in)
Net weight		23 kg (51 lbs)
Gross weight		29 kg (64 lbs)
Implementation according to EN 60601-1, EN 12021		Class I.
Classification acc. to MDD 93/42 EEC, 2007/47 EC		II b

a) The sound pressure level measured at a distance of 1.5 m using an A-weighting (filter).

Climatic conditions for storage and transport

Temperature -25°C to +55°C (-13°F to +131°F), 24 hrs +70°C (+158°F)

Relative air humidity 0% to 100% (with condensation)

Climatic conditions for operation

Temperature +15°C to +40°C (+59°F to +104°F)

Relative air humidity up to +70%

IPX0 Rating

FAD efficiency correction for differences in elevation

FAD correction table

Elevation [mamsl]	0 - 1500	1501 - 2500	2501 - 3500	3501 - 4500
FAD [l/min]	FAD x 1	FAD x 0.8	FAD x 0.71	FAD x 0.60

FAD efficiency refers to conditions at an elevation of 0 mamsl:

Temperature: 20°C

Atmospheric pressure: 101325 Pa

Relative humidity: 0%

OPERATION

4. OPERATION

INSTALLATION AND FIRST OPERATION

Do not use the compressor immediately after unpacking as it will not adjust to the ambient temperature.

Only qualified personnel can install the compressor and put it into operation for the first time. The installer shall train the operating personnel in the use and routine maintenance of the equipment. Installation and training of personnel should be acknowledged by the installer's signature on the installation certificate.

Before starting the compressor for the first time the stabilisation elements used to support the equipment during transport must be removed. Starting the compressor without removing the stabilisation elements can permanently damage the equipment!

The compressor does not contain a backup power supply.

It is forbidden to block the ventilation openings located in the upper part of the equipment!

If the compressor is equipped with a main source of air, the standby air source must be available.

ANY MODIFICATION OF THIS EQUIPMENT IS FORBIDDEN!


If this equipment is used nearby other instruments, the equipment must be observed in order to verify normal operations in the configuration it will be used.

Instruments may be affected electro-magnetically!

You may notice a "new product" odour when you first place the product into service (for a short period of time). This odour is temporary and does not impede the normal use of the product. Ensure the space is properly ventilated after installation.

Removal of transport stabilizers

Remove the four screws from the side of the equipment. Remove the cover and disconnect the ground wire.

Remove the two (2) M6 screws marked with the red warning washers and remove the stabilisation elements

Keep the screws and stabilisation elements in the event the compressor is transported at a later time.

Connect the grounding wire and reattach the cover.

Compressed air connection

The given appliance or respiratory apparatus should be connected to the **OUT** (1) quick coupler, as this is the outlet for compressed air.

The hose connecting the compressor to the respiration equipment must not pass through a cold environment i.e. placed on the ground. It should be as short as possible with no kinks (this may cause water to condense inside the hose).

The compressor must not be used for supplying air to the medical air central distribution pipeline system.

Electrical connection

The compressor comes with a plug containing an appropriate protective contact (ground.) Adhere to local electrical regulations. The voltage and frequency of the mains must comply with the specifications on the data label.

The electrical cord must not be stressed or have any tension exerted upon it, and must always be free.

OPERATION

- The socket must be accessible for safety reasons so that the equipment can be safely disconnected from the power supply in case of an emergency.
- The relevant current circuit must be protected.
- Connection of the ground connection (7) to other equipment must adhere to local electrical regulations.
- Fasten the electrical cord through the holder (6).

First operation

- Make sure that the stabilizing screws used during transport were removed.
- Check that the connection to the compressed air supply is correct.
- Check for proper connection to the main power supply.
- Switch on the pressure switch (5) to position "I".

The compressor will be operational once it is put into operation.

During equipment operation condensed water is released from the compressed air circuit through the automatic separator in the filter and drained to the evaporation tank.

Adjusting the output air pressure

The manufacturer has pre-set the output air pressure to 3 bar. If needed the output air pressure can be adjusted by turning the handle on the pressure regulator (2). The safety cover (15) must be removed before adjusting the output pressure. Push up the regulator knob and adjust the setting; once set, push it back down to secure its position. Reinstall the safety cover!

The maximum allowed output air pressure is 3.5 bar. Damage to respiration equipment may occur at higher outlet air pressures due to the presence of condensed water in the compressed air exiting the equipment.

OPERATION

THE EQUIPMENT MAY ONLY BE OPERATED BY THE TRAINED STAFF!

In case of emergency, switch the equipment off at the switch and pull out the main power plug.

Switching the compressor on

The compressor is switched on at the main power switch (1) by putting it in position "I".

Running the compressor

The compressor runs permanently once switched on. Condensed water removed from the compressed air evaporates from the evaporation tank back into the surrounding environment. The pressure gauge displays the outlet air pressure.

Cleaning the compressor

To clean the compressor, use a detergent that contains no abrasives, chemical solvents or other corrosive agents.

5. MAINTENANCE

REPAIRS AND SERVICE

Warranty and extended warranty repairs are to be completed by the manufacturer or a service provider authorized by the manufacturer.

The manufacturer reserves the right to modify the equipment in any way that will not alter the function or the operation of the equipment.

Only a qualified technician or the Customer Service Department of the manufacturer may perform repairs that go beyond routine maintenance. Use only spare parts and optional equipment approved by the manufacturer.

Prior to any maintenance or repairs, switch off the compressor and disconnect it from the mains (pull out the main power plug).

Air pump components (head, cylinder, pressure hose, etc.) are very hot during and shortly after compressor operation – do not touch these components!

Allow the equipment to cool down before maintenance, service or connecting/disconnecting the compressed air supply!

Automatic start: when pressure in the pressure tank decreases below the switch-on pressure, the compressor automatically switches on. The compressor automatically switches off when pressure in the air tank reaches the switch-off pressure.

Cover removal

- Remove the four screws from the side of the equipment.
- Remove the cover
- Disconnect the grounding conductor
- Reassemble using the opposite order

MAINTENANCE SCHEDULE

The operating entity is obliged to ensure that all tests of the equipment are carried out repeatedly at least once within every 24 months (EN 62353) or in intervals as specified by the applicable national legal regulations. A report must be prepared on the results of the tests (e.g.: according to EN 62353, Annex G), including the measurement methods used.

Maintenance	Page	Time interval	To be performed by
Replace filters in filter	13	Every 4000 hours	Qualified expert
Test the tightness of joints and inspect the equipment	16	After two years	Qualified expert
Replace piston complete with piston rings, O-rings and bearing	Service documentation	Every 8000 hours	Qualified expert
Replace suction filter	17	Every 8000 hours	Qualified expert
Perform "Repeated Test" according to EN 62353	15	1 x 2 years	Qualified expert

Checking tightness of joints and inspecting the equipment

Test for leakage:

- Disconnect the outlet hose from the quick coupling OUT (1); compressed air consumption will stop.
- Switch off the compressor
- Use soapy water to check for leaks in connections
- Tighten or re-seal joints as necessary.

Inspect the equipment:

- Check the status of the compressor motor for balanced operation or noise.
- Check the condition of the hangers above the pump.
- Check fan functionality.
- Check the working pressure values. Make sure the control knob on the pressure regulator (2) and the relief valve (11) are securely closed. The knob must be pulled up before setting the pressure or pushed down once the pressure has been set. Turn the pressure regulator (2) knob in the (+) direction until its maximum position is reached. Turn on the compressor and set the air consumption to around 10 L/min. Check the working pressure on the pressure gauge is at 5 bar; if this value needs to be changed, adjust the pressure using the knob on the relief valve (11). Turn

MAINTENANCE


this knob (2) to set the output pressure (normally 3 bar) at consumption of around 10 L/min.

- Check the condition of filter (12). With regular operation, condensate drains automatically. Verify this function by comparing the level of condensate in the tank to the scale line indicating the maximum level. If the level of condensate is above the scale line, replace the faulty parts.
- Examine the pump for:
 - Defects in the crankcase
 - Free movement of the crank shaft
- If necessary, replace any faulty parts.


Replacing filter elements

Prior to servicing the equipment, reduce the pressure to zero and disconnect the equipment from the main power supply.

- Unscrew the cover and pull it out using the enclosed wrench (16).
- Unscrew the filter holder.
- Replace the filter and screw the filter holder into place.
- Screw the cover back in place.

Spare filtration pad:

Plastic (5 µm) (12) AF20P-060S

Item 025200113-000

If the water release system becomes clogged, disassemble the release valve (f), clean the outlet, remove the float (b) with a lever mechanism (c, d, and e) from the tank, clean all components with a soap solution and then reassemble.

Changing suction filter elements

- Raise the clasps on the cover (9) and remove the cover.
- Replace the filter element. Filters that are only slightly dirty can be cleaned using compressed air. Filters can be washed.
- Reposition the cover.

Replacement parts:

Suction filter 03 POLYESTER, No. 025200126-000

Securing the equipment before transport

The equipment must be secured before any type of transport. Proceed using the opposite order of the procedure shown in Chapter 4 in the section for releasing the compressor.

SHUTDOWN

If the compressor is not going to be used for a long period of time, disconnected it from the main power supply and release the pressure.

EQUIPMENT DISPOSAL

- Disconnect the equipment from the main power supply.
- Release the pressure.
- Dispose of the equipment according to local regulations.
- Parts used in this product have no negative impact on the environment when disposed of properly.

TROUBLESHOOTING

6. TROUBLESHOOTING

Prior to servicing the equipment, reduce the pressure to zero and disconnect the equipment from the main power supply.

Only trained service personnel can perform the activities listed in the troubleshooting guide.

FAILURE	POTENTIAL CAUSE	SOLVING COMMON PROBLEMS
Compressor does not start	Problem with electrical power source	Main breaker is off
		Check voltage in socket
		Check fuses – replace faulty fuse
		Loosen terminal – tighten it
	Check the electrical cord - replace defective cord	
Compressor does not start	Damaged motor winding, damaged thermal protection	Replace the motor
	Capacitor failure	Replace capacitor
	Seized piston or other rotating part	Replace damaged parts
Low compressor performance	Incorrect working or output pressure settings	Adjust to stipulated pressure
	Air leak in compressed air distribution system	Check compressed air distribution system – seal loose joint
	Excessive air consumption	Do not exceed maximum flow (see Technical Data)
	Worn piston rings	Replace worn piston rings
	Plugged suction filter element	Replace contaminated element
	Plugged filter element	Replace contaminated element
Low compressor output pressure	Incorrect working or output pressure settings	Adjust to stipulated pressure
	Air leak in compressed air distribution system	Check compressed air distribution system – seal loose joint
	Excessive air consumption	Do not exceed maximum flow (see Technical Data)
	Worn piston assembly	Replace the piston assembly
	Plugged suction filter element	Replace contaminated element
	Plugged filter element	Replace contaminated element
Motor has stopped - thermal motor protection has tripped	Ventilation openings are covered	Uncover ventilation openings
	Cooling fans are not working	Replace the defective fan
	Equipment is hot, space is not ventilated	Secure suitable conditions
	Seized piston or other rotating part	Replace damaged parts
Compressor is noisy (knocking, metal noises)	Damaged piston bearing, piston rod, motor bearing	Replace damaged bearing
	Failed (cracked) hanger spring	Replace the damaged hanger
Water leaking from outlet	Hose from the equipment is too long or on the ground	Shorten the hose
	Malfunctioning filter float system	Clean or replace float system
	Incorrect working or output pressure settings	Adjust to stipulated pressure
	Excessive air consumption	Do not exceed maximum flow (see Technical Data)

7. SPARE PARTS


- Filtration elements AF20P-060S 025200113-000
- Fuse version 230V T6.3A 038100004-000
- Insertion DISS 1160-A 024000261-000
- suction filter 03 POLYESTER 025200126-000

ELECTRIC AND PNEUMATIC DIAGRAMS

8. ELECTRIC AND PNEUMATIC DIAGRAMS


WIRING DIAGRAM

1/N/PE ~ 230V, 50Hz
 ELEKTRIC OBJECT OF 1st. CAT.

EV1, EV2	Fan	Cb1	Capacitor
YV1	Solenoid valve	PH	Hour counter
M1	Elektric motor		

PNEUMATIC DIAGRAM

- | | |
|-------------------|-------------------------------------|
| 1. Suction filter | 6. Filter with condensate separator |
| 2. Compressor | 7. Safety valve |
| 3. Solenoid valve | 8. Pressure regulator |
| 4. Noise muffler | 9. Pressure gauge |
| 5. Cooler | 10. Compressed air outlet |

ELECTROMAGNETIC COMPATIBILITY DECLARATION

9. ELECTROMAGNETIC COMPATIBILITY DECLARATION

Medical equipment needs special precautions regarding electromagnetic compatibility (EMC) and needs to be installed and put into service according to the EMC information provided below.

Guidance and manufacturer's declaration - electromagnetic emissions		
Pursuant to IEC 60601-1-2:2014 - Medical electrical equipment - Part 1-2: General requirements for basic safety and essential performance - Collateral Standard: Electromagnetic disturbances		
The equipment is intended for use in the electromagnetic environment specified below. The customer or the user of the equipment should assure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	The equipment uses RF energy only for its internal functions. Therefore, the RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The equipment is suitable for use in all establishments, including domestic establishments and those directly connected to the public low voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	The equipment is not likely to cause any flicker, as the current flow is approx. constant after the start up.	

ELECTROMAGNETIC COMPATIBILITY DECLARATION

Guidance and manufacturer's declaration - electromagnetic immunity			
Pursuant to IEC 60601-1-2:2014 - Medical electrical equipment - Part 1-2: General requirements for basic safety and essential performance - Collateral Standard: Electromagnetic disturbances			
The equipment is intended for use in the electromagnetic environment specified below. The customer or the user of the equipment should assure that it is used in such an environment.			
Immunity test	IEC 60601-1-2 test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±8 kV contact ±15 kV air	±8 kV contact ±15 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	±2 kV for power supply lines ±1 kV for input/output lines	±2 kV 100 kHz repetition frequency Applied to mains connection	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±1 kV differential mode ±2 kV common mode	±1 kV L-N ±2 kV L-PE; N-PE Applied to mains connection	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruption, and voltage variations on power supply input lines IEC 60601-4-11	$U_T=0\%$, 0,5 cycle (at 0, 45, 90, 135, 180, 225, 270 and 315°) $U_T=0\%$, 1 cycle $U_T=70\%$ 25/30 cycles (at 0°) $U_T=0\%$, 250/300 cycles	$U_T \geq 95\%$, 0,5 cycle (at 0, 45, 90, 135, 180, 225, 270 and 315°) $U_T \geq 95\%$, 1 cycle $U_T=70\%$ (30% dipp U_T), 25(50Hz)/30(60Hz) cycles (at 0°) $U_T \geq 95\%$, 250(50Hz)/300(60Hz) cycle	Mains power quality should be that of a typical commercial or hospital environment. The device stops and restarts automatically at each dip. In this case does not occur unacceptable pressure drop.
Power frequency (50/60 Hz) IEC 61000-4-8	30 A/m	30 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment
NOTE: U_T is the A.C. mains voltage prior to application of the test level.			

ELECTROMAGNETIC COMPATIBILITY DECLARATION

Guidance and manufacturer's declaration - electromagnetic immunity			
Pursuant to IEC 60601-1-2:2014 - Medical electrical equipment - Part 1-2: General requirements for basic safety and essential performance - Collateral Standard: Electromagnetic disturbances			
The equipment is intended for use in the electromagnetic environment specified below. The customer or the user of the equipment should assure that it is used in such an environment.			
Immunity test	IEC 60601-1-2 test level	Compliance level	Electromagnetic environment - guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80MHz	3 Vrms	Portable and mobile RF communications equipment should be used no closer to any part of the equipment, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance $d=1,2\sqrt{P}$ $d=1,2\sqrt{P}$, 80 MHz to 800 MHz $d=2,3\sqrt{P}$, 800 MHz to 2,7 GHz
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2,7 GHz	3 V/m	where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m).
Proximity fields from RF wireless communications equipment IEC 61000-4-3	9 to 28 V/m 15 specific frequencies (380 to 5800 MHz)	9 to 28 V/m 15 specific frequencies (380 to 5800 MHz)	Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey ^a , should be less than the compliance level in each frequency range ^b . Interference may occur in the vicinity of equipment marked with the following symbol: 
NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies. NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			
^a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the equipment is used exceeds the applicable RF compliance level above, the equipment should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the equipment. ^b Over the frequency range 150 kHz to 80 kHz, field strengths should be less than 3 V/m.			

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