



DK50-10 S/M (MOBILE)



DK50 PLUS MOBILE



Accessories/ Zubehör/ Accessoires: **BOX** **FS** **REG** **ADK**

	DK50-10 S/ M (MOBILE)	
DRYER	✓	
BOX	✓	
(6,5-8bar)	50 l.min ³ / 7bar (50Hz)	60 l.min ³ / 7bar (60Hz)
S1-S3	S1 - 100%	
	≤ 58dB(A) (50Hz)	≤ 60dB(A) (60Hz)

	DK50 PLUS MOBILE	
DRYER	-	
BOX	✓	
(6-8bar)	70 l.min ³ / 6bar	
S1-S3	S1 - 100%	
	≤ 57dB(A)	