



DK50 4VR/50



DK50 2x4VR/110



Accessories/ Zubehör/ Accessoires: **BOX** **MD** **FS** **REG** **AOK**

	DK50 4VR/50	DK50 4VR/50/M	DK50 4VR/50S	DK50 4VR/50S/M
DRYER	-	✓	-	✓
BOX	-	-	✓	✓
(6-8bar)	250 l.min ⁻¹ / 6bar	200 l.min ⁻¹ / 6bar	250 l.min ⁻¹ / 6bar	200 l.min ⁻¹ / 6bar
(8-10bar)	210 l.min ⁻¹ / 8bar	140 l.min ⁻¹ / 8bar	210 l.min ⁻¹ / 8bar	140 l.min ⁻¹ / 8bar
S1-S3	S1 - 100%	S1 - 100%	S1 - 100%	S1 - 100%
	≤ 77 dB(A)	≤ 77 dB(A)	≤ 58 dB(A)	≤ 58 dB(A)

Accessories/ Zubehör/ Accessoires: **BOX** **MD** **FS** **REG**

	DK50 2x4VR/110	DK50 2x4VR/110/M	DK50 2x4VR/110S	DK50 2x4VR/110S/M
DRYER	-	✓	-	✓
BOX	-	-	✓	✓
(6-8bar)	520 l.min ⁻¹ / 6bar	430 l.min ⁻¹ / 6bar	520 l.min ⁻¹ / 6bar	430 l.min ⁻¹ / 6bar
(8-10bar)	410 l.min ⁻¹ / 8bar	330 l.min ⁻¹ / 8bar	410 l.min ⁻¹ / 8bar	330 l.min ⁻¹ / 8bar
S1-S3	S1 - 100%	S1 - 100%	S1 - 100%	S1 - 100%
	≤ 77	≤ 77	≤ 58	≤ 58