



## DK50 PLUS



## DK50 2V



Accessories/ Zubehör/ Accessoires: **BOX MD FS REG AOK**

	DK50 PLUS	DK50 PLUS/M	DK50 PLUS S	DK50 PLUS S/M
DRYER	-	✓	-	✓
BOX	-	-	✓	✓
(6-8bar)	70 l.min <sup>-1</sup> / 6bar (50Hz) 80 l.min <sup>-1</sup> / 6bar (60Hz)	60 l.min <sup>-1</sup> / 6bar (50Hz) 70 l.min <sup>-1</sup> / 6bar (60Hz)	70 l.min <sup>-1</sup> / 6bar (50Hz) 80 l.min <sup>-1</sup> / 6bar (60Hz)	60 l.min <sup>-1</sup> / 6bar (50Hz) 70 l.min <sup>-1</sup> / 6bar (60Hz)
(8-10bar)	60 l.min <sup>-1</sup> / 8bar (50Hz) 70 l.min <sup>-1</sup> / 8bar (60Hz)	50 l.min <sup>-1</sup> / 8bar (50Hz) 60 l.min <sup>-1</sup> / 8bar (60Hz)	60 l.min <sup>-1</sup> / 8bar (50Hz) 70 l.min <sup>-1</sup> / 8bar (60Hz)	50 l.min <sup>-1</sup> / 8bar (50Hz) 60 l.min <sup>-1</sup> / 8bar (60Hz)
S1-S3	S1 - 100%	S1 - 100%	S1 - 100%	S1 - 100%
	≤ 66 dB(A)	≤ 66 dB(A)	≤ 47 dB(A)	≤ 47 dB(A)

Accessories/ Zubehör/ Accessoires: **BOX MD FS REG AOK**

	DK50 2V	DK50 2V/M	DK50 2VS	DK50 2VS/M
DRYER	-	✓	-	✓
BOX	-	-	✓	✓
(6-8bar)	135 l.min <sup>-1</sup> / 6bar	104 l.min <sup>-1</sup> / 6bar	135 l.min <sup>-1</sup> / 6bar	104 l.min <sup>-1</sup> / 6bar
(8-10bar)	110 l.min <sup>-1</sup> / 8bar	87 l.min <sup>-1</sup> / 8bar	110 l.min <sup>-1</sup> / 8bar	87 l.min <sup>-1</sup> / 8bar
S1-S3	S1 - 100%	S1 - 100%	S1 - 100%	S1 - 100%
	≤ 67 dB(A)	≤ 69 dB(A)	≤ 51,5 dB(A)	≤ 55 dB(A)