



DK50 SIMPLE



	DK50 simple
DRYER	-
BOX	✓
(4,5-6bar)	50 l.min ⁻¹ / 5bar
S1-S3	S1 - 100%
	≤ 53 dB(A)

DK50 B



Accessories/ Zubehör/ Accessoires: **BOX** **FS** **REG** **AOK**

	DK50 B	DK50 BS
DRYER	-	-
BOX	-	✓
(9,5-11,5bar)	50 l.min ⁻¹ / 10bar	
S1-S3	S3 - 50%	S3 - 50%
	≤ 65 dB(A)	≤ 45 dB(A)