

# DECLARATION OF MANUFACTURER

No.: PV-1/2010/RQ

**Manufacturer:** Ekom spol. s r.o.

**Address:** Priemysel'na 5031/18  
921 01 Piešťany  
Slovak Republic  
ID No. 31416519

**Product:**

Product Name:  
Type:

**Dental Compressor**

DK50 B, DK50 BS, DK50 I, DK50 H, DK50C, DK50 K, DK50 P, DK50 M, DK50 Z, DK50 Z/K, DK50 S, DK50 S/K, DK50-10 I, DK50-10 Z, DK50-10 Z/K, DK50-10 S, DK50-10 S/K, DK50-10 Z/M, DK50-10 S/M, DK50 PLUS, DK50 PLUS/K, DK50 PLUS/M, DK50 PLUS S, DK50 PLUS S/K, DK50 PLUS S/M, DK50 2V, DK50 2V/K, DK50 2V/M, DK50 2VS, DK50 2VS/K, DK50 2VS/M, DK50 2V/110, DK50 2V/110/K, DK50 2V/110/M, DK50 2V/110S, DK50 2V/110S/K, DK50 2V/110S/M, DK50 2x2V/110, DK50 2x2V/110/K, DK50 2x2V/110/M, DK50 2x2V/110S, DK50 2x2V/110S/K, DK50 2x2V/110S/M, DK50 4x2V/M, DK50 4x2V, DK50 6x2V/M, DK50 6x2V, DK50 4x2VT/M, DK50 6x2VT/M, DK50 9x2VT/M, DK50 12x2VT/M, DK50 15x2VT/M, DK50 4x2VT, DK50 6x2VT, DK50 9x2VT, DK50 12x2VT, DK50 15x2VT, DO 2.1, DO 2.1T, DO 2.1Z, DO 2.1TZ, DO 2.1-10, DO 2.1-10T, DO 2.1-10Z, DO 2.1-10TZ, DO-M, DUO, DUO 2V, DUO/M, DUO 2V/M, DUO T, DUO T/M, DUO 2VT, DUO 2VT/M, DUO 2, DUO 2/M, DUO 2T, DUO 2T/M

The manufacturer declares that the manometer used in the above compressor types is of indicative character and serves only to inform the operators of the presence of pressure in the device.

The manometer is not intended for the measurement of pressure in a pneumatic system. It is necessary to use an independent calibrated manometer with the required accuracy for the purposes of setting the elements in the pneumatic system.

Piešťany 23.4.2010

 - 26 -  
 spol. s r.o.  
 Priemysel'na 5031/18  
 921 01 Piešťany  
 IČO: 31416519 IČ DPH: S102029390366

  
 .....  
 Ing. Vladimír Batora  
 Director

