

# REGULAR MAINTENANCE AND SERVICE REPAIR RECORD



Information about maintenance or repair – please tick performed activities

Product serial number:	
Type of product:	
Hour counter value ( if available) OR Number of years of operation (since installation):	

## A. Regular maintenance – please tick performed activities

Recommended service intervals (following specified working environment conditions in user manual)	Without hour counter			With hour counter
	by every inspection	every year	every 2 years	every 2000 hours
Equipment cleaning (cooler, fan) <sup>1</sup>	<input type="checkbox"/>	<input type="checkbox"/>	-	<input type="checkbox"/>
Leakage test (by soapwater ) <sup>1</sup>	<input type="checkbox"/>	<input type="checkbox"/>	-	<input type="checkbox"/>
Inspection of el. connections <sup>1</sup>	<input type="checkbox"/>	<input type="checkbox"/>	-	<input type="checkbox"/>
Safety valve check <sup>1</sup>	-	<input type="checkbox"/>	-	<input type="checkbox"/>
Replace inlet filter <sup>1</sup>	-	-	<input type="checkbox"/>	<input type="checkbox"/>
Replace inlet pre-filter (2V, 4VR compressors) <sup>1</sup>	-	-	<input type="checkbox"/>	<input type="checkbox"/>
Replace membrane dryer filter and micro filter elements (for compressors equipped with dryer) <sup>1</sup>	-	<input type="checkbox"/>	-	<input type="checkbox"/>

## B. Repair service (please list replaced parts and their amount)

Name and article number of the spare part	Pcs.

## C. General evaluation of the product after service (please describe):

Date of maintenance/ repair:		Name of service company:	
Service engineer name and surname:			

<sup>1</sup>/for detailed reference of particular action please check user manual

<sup>2</sup>/other language versions of this report are available on Ekom web page in Service-Trade zone.

<sup>3</sup>/if more than in booklet provided quantity of reports are required, please make a copy or download it from Ekom web page - section Forms in Service-Trade zone