

DO 3 DENTAL SUCTION SYSTEM



Compatible with all wet or semiwet dental units

- ✓ Solution for 1 up to 4 dental units
- ✓ Uninterrupted operation during aspiration
- ✓ Quietest solution on the market thanks to noise reduction cabinet
- ✓ Touch screen control panel
- ✓ Service warnings and operation logs
- ✓ Build in function test of key components for easy troubleshooting
- ✓ Optional outlet HEPA filter



Touch screen control panel



Easy maintenance



Technical specification	DO 3			
	BASIC	STANDARD	ADVANCED	
Model	BASIC	STANDARD	ADVANCED	
Performance	300	600	900	lit.min-1
Nominal voltage, frequency	230, 50/ 60	230, 50/ 60	230, 50/ 60	V, Hz
Number of connected dental units	1	1 - 2	3 - 4	
Coefficient of use	100	100	75	%
Maximum vacuum	-130	-160	-200	mbar
Noise level with/ without cabinet	54/ 60	55/ 62	-/ 64	dB(A)
Operating mode	intermittent	continuous	continuous	
Dimensions without the amalgam separator - without cabinet / with cabinet				
Height	520 / 545	520 / 545	590	mm
Depth	420 / 445	420 / 445	540	mm
Width	360 / 405	360 / 405	360	mm
Weight	21 / 32	25 / 35	41	kg

Model	BASIC 300 l/min 1 dental chair	STANDARD 600 l/min 1-2 dental chairs	ADVANCED 900 l/min 3-4 dental chairs	SOUNDPROOFING CABINET	CONNECTION FOR AS*
DO 3 Z BASIC	✓				
DO 3 BASIC	✓			✓	
DO 3 ZA BASIC	✓				✓
DO 3 A BASIC	✓			✓	✓
DO 3 Z STANDARD		✓			
DO 3 STANDARD		✓		✓	
DO 3 ZA STANDARD		✓			✓
DO 3 A STANDARD		✓		✓	✓
DO 3 Z ADVANCED			✓		
DO 3 ZA ADVANCED			✓		✓

*AS - sedimentation amalgam separator



EKOM spol. s r.o., Priemysel'na 5031/18, 921 01 Piešťany, Slovak Republic
T +421 33 7967 205, sales@ekom.sk, www.ekom.sk

